



CELOX™

Control Bleeding. Save Lives.

CELOX-A controls bleeding of penetrating wounds.

Celox-A™ is a new unique applicator delivery system design to help control life-threatening bleeding from small penetrating wounds.

Hemostatic agents are not easily applied through small entry wounds such as wounds from shrapnel, blades or penetrating trauma. Celox-A™ allows the hemostatic agent to be applied to a deep bleeding site using its applicator and plunger system.

Each package contains an applicator system pre-filled with 6g of Celox™ .



FEATURES:

- Helps save lives by controlling bleeding within minutes.
- Works in hypothermic conditions.
- Clots blood containing Coumadin / Warfarin.
- Clots blood containing Heparin.
- No heat—No Burning.
- Sterile waterproof packaging.
- No special requirements.
- Easy to remove from the wound.

The Celox-A™ applicator system is the subject of a PCT patent application.

QUOTES From the field.

- CELOX is without a doubt the most exciting hemostatic dressing I have had the opportunity to work with. It is extremely easy to use, almost foolproof. I believe every soldier, not just every medic should carry CELOX on their kit. (SSG Medic).
- I finally got to try CELOX, and was very impressed. Had a casualty in a controlled environment with deep femoral bleed. I was very impressed with your product . 68W (Combat Medic) Instructor .
- The CELOX product not only surprised me but gave me hope that this may make it to our troops and since I noticed no heat generation by the CELOX product and the fact that even with improper application it works very well, that if applied correctly with a very serious injury it will perform as expected without the side effects....." Fire Dept. - TX.
-I had the opportunity to use the CELOX numerous times on some serious wounds and was very impressed. I found it to be the easiest to use and the easiest to train people to use. It also allowed for the correction of improper placement after the initial application. (Retired 18-D)

