



CELOX™

Be prepared... Celox Saves Lives.

CELOX-A gains prestigious class 3 CE mark

CELOX-A™ FIELD REPORT.


The Celox applicator that I used was very easy. I actually primed the end of the barrel a little to make sure that it was not clogged. I inserted it about 2 inches into the wound and inserted all of its contents without any problems in about 1 second.

The wound was a blast injury with some major distal bleeding. The laceration was about 7 inches in length and full thickness to the bone with major cavitation.

One of the distal vessels in the leg was bleeding very badly and when I saw him he had on a tourniquet. The Celox was introduced with a bandage and the tourniquet was removed. There was no further bleeding from the wound.

This was my first time using the new applicator and as I expected, it worked amazingly. The tube filled the wound deep and directly where I needed the Celox. Of course the Celox did an amazing job getting the bleeding under control and made the patient more stable, manageable and was a life saving intervention.'

US Special Forces Medic



MedTrade, Crewe, UK is proud to announce that Celox-A, the unique life saving product, has obtained EU Class CE Mark approval .

Celox-A™ is a unique applicator delivery system designed to stop life-threatening bleeding from a range of wounds including difficult to treat penetrating trauma caused by shrapnel, knives or bullet entry.

Penetrating trauma often results in serious bleeding from vessels deep in the body where the initial emergency treatment has to be given through a small entry wound. Such wounds are often fatal. Celox-A™ allows the Celox™ haemostat to be applied to a deep bleeding site using its applicator system. It makes such treatment easy, safe and very effective.

Since its introduction to the US market in 2006 Celox's™ ability to rapidly clot blood has made it the haemostat of choice for many US military, police, fire and ambulance medics and now also for workplace first aid. Celox™ has been demonstrated to rapidly and reliably stop major bleeding. It works independently of body temperature and is able to clot blood in the presence of anti-coagulants such as heparin and Coumadin. CE approval provides European military, police, emergency and paramedic teams with an invaluable tool to help save lives.

Craig Hardy, CEO of MedTrade, said " Celox-A™ is an outstanding medical device. It has been a true pleasure to hear of the lives it has already saved in Afghanistan. We are now proud that Celox-A™ has obtained a class 3 CE Mark. Class 3 is the highest grade of medical device and approval is only given after a very thorough review of a product's safety and performance by an independent expert authority. A class 3 CE mark is a real mark of quality."

For more details www.CeloxMedical.com or Celox@MedTrade.co.uk.

Celox™ & Celox-A™ were developed by MedTrade Products Ltd, Crewe, UK, (www.MedTrade.co.uk). There are opportunities for overseas distribution.

Celox-A will be at the Medica, Dusseldorf from the 19 to 22 Nov 08. Visit us at Hall 16 F18-3.

Email Celox@medtrade.co.uk for an appointment.